

ABT RINSING



ABT RINSING SOLUTION is used for rinsing embryo collection filters after the collection procedure is completed. ABT's formulation is the same as our ABT Complete Flush but excludes PVA to help avoid foaming during use making embryo searching more efficient by not having to wait for the bubbles to clear. Using the same basic formula from flushing to rinsing means that embryos are kept in the same base solution as they were collected in.

Using ultra-purified water as a base, this solution is filtered to 0.2uM so there is no need for additional filtering before use. ABT products are all manufactured in a cleanroom environment and go through strict quality control practices ensuring consistency and quality from lot to lot.

- Excludes PVA to avoid foaming
- Contains Gentamycin as an antibiotic
- Animal origin free & 100% traceable components
- Tight quality control standards equal consistent results
- Durable EVA packaging with a spike port and non-latex needle septum removable cap
- Made in the USA with two decades of experience

PRODUCT # SIZE
ABT100 **1 LITER**

Recommended Storage: 15 - 30 °C, however, it will not be harmed by refrigeration or freezing.

Shelf-life: 24 months from date of manufacture.

ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H
PULLMAN, WA 99163 USA
Phone: 509-592-8144
www.abt360llc.com
alyssa@abt360llc.com
facebook.com/ABT360
instagram.com/ABT_360