



1615 NE Eastgate Blvd, Sec H
Pullman, WA 99163 USA
Phone 509.592.8144
www.abt360llc.com

ABT COMPLETE FLUSH

ABT014, 1 LITER
ABT015, 2 LITER
ABT018, 4 LITER

This solution is “ready to use” and free of any components of animal origin. Heavy metal chelators have been added and Polyvinyl alcohol (PVA) is the surfactant.

Recommended storage: 15 - 30 °C; however, it will not be harmed by refrigeration or freezing.

ABT RINSING

ABT100, 1 LITER

Recommended Use: To help ensure all embryos have been retrieved from the filter, use a 20ml syringe filled with ABT Rinsing Solution to wash the collection filter membrane. Hold the filter over a clean searching dish and use the syringe to spray solution over the filter membrane.

Recommended storage: 15 - 30 °C; however, it will not be harmed by refrigeration or freezing.

Disclaimer

The recommended procedures in this insert are not considered a complete protocol and will vary according to user. ABT media are perishable products. Storage and handling will influence their rate of deterioration. ABT 360, LLC does not warrant any viability or conception resulting from the use of an ABT product. For further information contact us at (509) 592-8144 or visit our website at www.abt360llc.com.