



1615 NE Eastgate Blvd, Sec H
Pullman, WA 99163 USA
Phone 509.592.8144
www.abt360llc.com

ABT TL HEPES

ABT070, 500ML

ABT TL Hepes is a medium used for rinsing oocytes in certain steps of IVF protocols. TL Hepes does not contain BSA. This product is packaged in sterile 500ml high-density PETG bottles. Do not use if product is not translucent or the tamper proof cap has been compromised.

Recommended Storage: 2 - 30°C. TL Hepes can be stored at room temperature however it will not be harmed by refrigeration. Do not freeze.

Shelf Life: 18 months

Disclaimer

The recommended procedures in this insert are not considered a complete protocol and will vary according to user. ABT media are perishable products. Storage and handling will influence their rate of deterioration. ABT 360, LLC does not warrant any viability or conception resulting from the use of an ABT product. For further information contact us at (509) 592-8144 or visit our website at www.abt360llc.com.